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Central Sterilization (Supply) Room :CSR ←

Central Sterile (Sterilization) Supply (Services) Department :CSSD ←

Central Sterilization (Service) Department :CSD ←

Theatre Sterile Supply Unit :TSSU ←

Sterile Processing Department (Distribution) :SPD ←

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CENTRAL STERILE:
THE HEART OF
THE HOSPITAL

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**NO ENTRANCE
AUTHORIZED
PERSONNEL ONLY**

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

SuperNova[®]
InstruCrema Lubricant Concentrate

Surgical Instrument Lubricant

Directions for Use: Designed for the final rinse cycle of automated washers. May be applied manually to the working parts of instrumentation or dilute for soaking. Use lubricant after a thorough cleaning and rinsing. Allow to dry; then proceed to autoclave.

Recommendation: May be added to the final, high purity rinse cycle of automated washes at the dilution rate below.

Temperature: May be used at room temperature for soaking and direct application.

Use Dilution: Soaking 1:10 proportion. Automated washer 1/10 fl. oz/gal. to 1/4 fl. oz/gal. (0.78 mL/L to 1.96 mL/L). Dilution based on water quality and machine setting.

Specially formulated aqueous emulsion designed to ease stiffness, prevent rust and extend the life of surgical instruments.

- Maintains smooth working action.
- Rapidly breaks down bioburden and organic soil
- Natural odor, milky color
- Water dispersible

- Validated pH neutral
- Environmentally preferred, bio-degradable
- US EPA/SAFER CHOICE recognition for sustainability and safer chemistry.

Order Number: CSNL01
Net Contents: 2.5 Gallon (9.46 liters)

SAFER CHOICE
Meets EPA's Safer Choice Criteria

Case Medical, Inc.
24 Empire Blvd., 5, Hackensack, NJ 07601
(888)227-CASE www.casemed.com

Storage conditions: Store at room temperature away from heat source. Keep container closed when not in use.
Compatibility: Non-corrosive, pH neutral and lower toxicity. Compatible with stainless steel, aluminum, titanium, carbon, copper, and brass.
Caution: Ingredients may cause severe skin burns and eye damage. Harmful if swallowed.
Hazard class not classified: Corrosive liquid, skin sensitizer >20%, Sodium Hydroxide <1%, Potassium Hydroxide <1%.

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

SEQUENCE FOR PUTTING ON PERSONAL PROTECTIVE EQUIPMENT (PPE)

The type of PPE used will vary based on the level of precautions required, such as standard and contact, droplet or airborne infection isolation precautions. The procedure for putting on and removing PPE should be tailored to the specific type of PPE.

1. GOWN

- Fully cover torso from neck to knees, arms to end of wrists, and wrap around the back
- Fasten in back of neck and waist

2. MASK OR RESPIRATOR

- Secure ties or elastic bands at middle of head and neck
- Fit flexible band to nose bridge
- Fit snug to face and below chin
- Fit-check respirator

3. GOGGLES OR FACE SHIELD

- Place over face and eyes and adjust to fit

4. GLOVES

- Extend to cover wrist of isolation gown

USE SAFE WORK PRACTICES TO PROTECT YOURSELF AND LIMIT THE SPREAD OF CONTAMINATION

- Keep hands away from face
- Limit surfaces touched
- Change gloves when torn or heavily contaminated
- Perform hand hygiene


HOW TO SAFELY REMOVE PERSONAL PROTECTIVE EQUIPMENT (PPE) EXAMPLE 1

There are a variety of ways to safely remove PPE without contaminating your clothing, skin, or mucous membranes with potentially infectious materials. Here is one example. **Remove all PPE before exiting the patient room** except a respirator, if worn. Remove the respirator **after** leaving the patient room and closing the door. Remove PPE in the following sequence:

1. GLOVES

- Outside of gloves are contaminated!
- If your hands get contaminated during glove removal, immediately wash your hands or use an alcohol-based hand sanitizer
- Using a gloved hand, grasp the palm area of the other gloved hand and peel off first glove
- Hold removed glove in gloved hand
- Slide fingers of ungloved hand under remaining glove at wrist and peel off second glove over first glove
- Discard gloves in a waste container

2. GOGGLES OR FACE SHIELD

- Outside of goggles or face shield are contaminated!
- If your hands get contaminated during goggle or face shield removal, immediately wash your hands or use an alcohol-based hand sanitizer
- Remove goggles or face shield from the back by lifting head band or ear pieces
- If the item is reusable, place in designated receptacle for reprocessing. Otherwise, discard in a waste container

3. GOWN

- Gown front and sleeves are contaminated!
- If your hands get contaminated during gown removal, immediately wash your hands or use an alcohol-based hand sanitizer
- Unfasten gown ties, taking care that sleeves don't contact your body when reaching for ties
- Pull gown away from neck and shoulders, touching inside of gown only
- Turn gown inside out
- Fold or roll into a bundle and discard in a waste container


4. MASK OR RESPIRATOR

- Front of mask/respirator is contaminated — DO NOT TOUCH!
- If your hands get contaminated during mask/respirator removal, immediately wash your hands or use an alcohol-based hand sanitizer
- Grasp bottom ties or elastics of the mask/respirator, then the ones at the top, and remove without touching the front
- Discard in a waste container

5. WASH HANDS OR USE AN ALCOHOL-BASED HAND SANITIZER IMMEDIATELY AFTER REMOVING ALL PPE

PERFORM HAND HYGIENE BETWEEN STEPS IF HANDS
BECOME CONTAMINATED AND IMMEDIATELY AFTER
REMOVING ALL PPE


HOW TO SAFELY REMOVE PERSONAL PROTECTIVE EQUIPMENT (PPE) EXAMPLE 2

Here is another way to safely remove PPE without contaminating your clothing, skin, or mucous membranes with potentially infectious materials. **Remove all PPE before exiting the patient room** except a respirator, if worn. Remove the respirator **after** leaving the patient room and closing the door. Remove PPE in the following sequence:

1. GOWN AND GLOVES

- Gown front and sleeves and the outside of gloves are contaminated!
- If your hands get contaminated during gown or glove removal, immediately wash your hands or use an alcohol-based hand sanitizer
- Grasp the gown in the front and pull away from your body so that the ties break, touching outside of gown only with gloved hands
- While removing the gown, fold or roll the gown inside-out into a bundle
- As you are removing the gown, peel off your gloves at the same time, only touching the inside of the gloves and gown with your bare hands. Place the gown and gloves into a waste container

2. GOGGLES OR FACE SHIELD

- Outside of goggles or face shield are contaminated!
- If your hands get contaminated during goggle or face shield removal, immediately wash your hands or use an alcohol-based hand sanitizer
- Remove goggles or face shield from the back by lifting head band and without touching the front of the goggles or face shield
- If the item is reusable, place in designated receptacle for reprocessing. Otherwise, discard in a waste container


3. MASK OR RESPIRATOR

- Front of mask/respirator is contaminated — DO NOT TOUCH!
- If your hands get contaminated during mask/respirator removal, immediately wash your hands or use an alcohol-based hand sanitizer
- Grasp bottom ties or elastics of the mask/respirator, then the ones at the top, and remove without touching the front
- Discard in a waste container

4. WASH HANDS OR USE AN ALCOHOL-BASED HAND SANITIZER IMMEDIATELY AFTER REMOVING ALL PPE

**PERFORM HAND HYGIENE BETWEEN STEPS IF HANDS
BECOME CONTAMINATED AND IMMEDIATELY AFTER
REMOVING ALL PPE**

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
FIBEROPTIC CABLE CARE

Recommended:

- Coil fiberoptic cable in a six inch diameter groove
- When cleaning storage material (boxes), use a soft brush with lukewarm water and a mild non-toxic soap solution
- Rinse after cleaning with distilled water
- Keep cables away from sharp instruments, edges, etc. Any case or cover can make the cable smooth. It should be taken out of service immediately
- For maximum comfort and maneuverability, use the green clips to secure the headlamp cable to the surgeon's gown

Not Recommended:

- Do not pull, stretch or kink fiberoptic cables
- Do not clean with strong soaps or solvents
- Do not use on reusable endoscopes
- Do not allow cable to strike floor or other hard surface
- Do not use on steel cables in a slightly curved diameter
- Do not use the cable's original design of lubrication. Do note the bundle size has been marked on the instrument aperture for maximum light output. NOTE: A larger aperture cable size will not increase the light output of a smaller aperture instrument.

CLEANING AND STERILIZATION

The fiberoptic cable should be cleaned of tissue and debris prior to sterilization in a neutral pH based detergent solution. Cable may be washed in disinfecting solution. The time for steam sterilization varies depending on the type of steam sterilant and the temperature at which it is run and whether the items are wrapped in paper or single material. Recommendations are as follows:

| Type of Steam Sterilant | Temperature | Time |
|--|---------------|------------|
| Gaseous Ethylene Oxide (wrapped) | 32°C (90°F) | 36 Minutes |
| | 52°C (125°F) | 22 Minutes |
| Steam-Flask Pressure Pulse (unwrapped) | 121°C (250°F) | 28 Minutes |
| Gaseous Ethylene Oxide (wrapped) | 32°C (90°F) | 18 Minutes |
| Bromine peracetic (wrapped) | 132°C (270°F) | 3 Minutes |

The following sterilizer types are recommended for use with cables:

- Ethylene Oxide (ETO)
- Steam

To avoid contact with sharp or pointed objects during purification and storage, use a sterilization bin. Make sure your sterilant is working properly and that specific instructions and exposure periods are followed. Any small holes that may appear in the sleeve during steam sterilizing will not harm the cable and will not affect its design.

After sterilizing, allow the cable to cool slowly to room temperature. CAUTION: Exposure to cold water will cause fiber breakage and a reduction of light output.

Manufactured for:

Teleflex
MEDICAL

Research Triangle Park, NC 27709 USA
USA (866) 246-6990 • International (919) 544-6000
220002432

Rx ONLY

CE
80299
T273-00000-1M
1491 10/2006 14912 2011 M

Keep Advancing and Enjoy Your Life as a Champion

IFU

Appendix A: Special Cleaning Considerations Table

SPECIAL CLEANING CONSIDERATIONS — HANDPIECES

| Component | Moving Parts | Critical Areas | Brush Sizes |
|--------------------------------|--------------|--|---------------------------|
| RemB Electric Micro Drill | Yes | Lock sleeve, ring, offset pin | 6.35 mm |
| Oscillating Saws | Yes | Blade mount, blade mount collar | 1 mm, 10 mm |
| Sagittal Saws | Yes | Blade mount head | 1 mm, 10 mm |
| Reciprocating Saws | Yes | Blade receptacle, button | 3.33 mm |
| CORE Universal Driver/U Driver | Yes | Lumen | 10 mm, 12.7 mm |
| Hummer/ESSx Handpieces | Yes | Suction lumen, cutter receptacle | 3.33 mm, 15 mm |
| Impaction Drill | Yes | Dead-ended lumen | 3 mm |
| RemB Electric Universal Driver | Yes | Attachment connector, lumen, triggers, safety bar, back cap crevices | 6.35 mm |
| RemB Electric Wiredriver | Yes | Lumen, lever, trigger, handle trigger recess, back cap crevices | 2.5 mm, tooth-brush style |
| CORE Micro Driver | Yes | Lumen | 6.35 mm |
| TPS Micro Driver | Yes | Lumen | 4 mm, 6 mm, 8 mm |

IFU

DePuy Synthes

Equipment: Ultrasonic cleaner, washer/disinfector, various sized soft-bristled brushes, lint-free cloths, syringes, pipettes and/or water jet, enzymatic cleaner or detergent solution

Pre-clean method (Pre-clean method must be performed prior to mechanical washer method listed below.)

1 MIN
+ 2 MINS

+ 1 MIN

1. Rinse soiled device under running cold tap water for a **minimum of one minute**. Remove gross soil using a soft-bristled brush or soft, lint-free cloth.
2. Manually clean device for a **minimum of two minutes** in a freshly prepared enzymatic cleaner or detergent solution. Follow the enzymatic cleaner or detergent manufacturer's instructions for the correct dilution, temperature, water quality and exposure time. Use a soft-bristled brush to remove soil and debris. Actuate joints, handles and other movable device features to expose all areas to detergent solution, if applicable. Clean device under water to prevent aerosolization of contaminants, *Note: fresh solution is a newly-made, clean solution.*
3. Rinse device using cold to lukewarm running tap water for a **minimum of one minute**. Use a syringe, pipette or water jet to flush lumens and channels. Actuate joints, handles and other moveable device features in order to rinse thoroughly under running water, if applicable.
4. Prepare a fresh detergent solution using enzymatic cleaner or detergent. Follow the enzymatic cleaner or detergent manufacturer's instructions for the correct dilution, temperature, water quality and exposure time. *Note: fresh solution is a newly-made, clean solution.*
5. Clean Synthes devices ultrasonically for a **minimum of 15 minutes**, using a minimum frequency of 40 KHz.
6. Rinse device using DI or PURW water for a **minimum of two minutes**. Use a syringe, pipette or water jet to flush lumens and channels. DI or PURW water must be used for final rinse.
7. Visually inspect device. Repeat steps 2–7 until no visible soil remains on device.

Cleaning – **15 MINS**
Automated/Mechanical washer Method

+ 2 MINS

**= ABSOLUTE
MINIMUM -
21 TOTAL MINS
OF PRECLEANING**

Mechanical Washer process: (Pre-cleaning steps 1–7 should occur prior to this step.) *Note: The washer/disinfector should fulfill requirements specified in ISO 15883. Use MIS injector unit to process lumens and cannulations.*

8. Process device using the following cycle parameters:

INSTRUCTIONS

Use

Before initial use: Remove the instruments from their packaging and remove the protective sheaths and protective wrappers. Before initial use, please read the following additional instructions.

After use: Disassemble the instruments. If they come into contact with corrosive medications or cleaning agents, wash them immediately in water. If there are encrustations on the instruments after cleaning, further chemical treatment is required (e.g. 3% H₂O₂ for 5 minutes).

Dry disposal: Dry disposal can be used for the instruments. Process the instruments as quickly as possible after disposal.

Preparation for decontamination

Wet disposal: The instrument must be completely immersed in a non-corrosive and non-fixing disinfectant bath, whereby the cleaning and disinfectant solutions must be replaced when they are visibly contaminated or at least daily.

Before cleaning, decontamination agents must be removed with water.

Cleaning and disinfection, manual

Mechanical cleaning and disinfection is to be given preference over manual processing.

Open the instruments in accordance with the instructions. Cleaning and disinfection solutions must be replaced if they are visibly contaminated or at least daily.

Ultrasound bath: After disinfection the instruments are placed in sieves so that no shadows are formed. Add cleaning agent to the water and adjust the temperature of the solution to comply with the instructions of the cleaning agent manufacturer. The cleaning in the ultrasound bath (35-40 kHz) takes 5 minutes.

Manual cleaning: Use lint-free cloths, synthetic brushes or water-jet guns for manual cleaning with a cleaning agent and disinfectant. Brush structured surfaces while they are immersed in the cleaning agent. Rinse with completely desalinated water and dry the instruments after the rinsing process. Moveable parts must be moved at least 20 times during manual processing. 

Cleaning and disinfection, automatic

Cleaning and disinfection unit: Open the cleaning cart and place the instruments gently and securely inside. Do not overload the cleaning cart. The preliminary rinse is followed by chemical cleaning at 40 to 60 °C. Residues from the cleaning process must be reliably removed in the final rinsing phase. Prevent

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Surgical Instrument Data

INVENTORY DATA

- Quantity of Assets
- Locate in Instruments in Trays
- Spare replacement inventory

STERILIZATION
4min 270F 40min Dry

WASH
131°F (55°C) with a 0.8% enzymatic detergent

| | |
|--------|-------|
| Usage | 350 |
| Price | \$275 |
| Repair | \$35 |
| TCO | \$385 |
| SKU | 120 |

PRE-CLEANING INSTRUCTIONS

- Soak in lukewarm tap water for 5 min.
- Flush with a 20ml syringe.
- Brush under lukewarm tap water with a soft brush until all visible residues are removed.
- Sonicate for 10 min, ultrasonic bath at 104°F (40°C) with 0.8% enzymatic cleaner.

TISCHLER BIOPSY PUNCH FORCEPS
CROCODILE ACTION 210MM EFFECTIVE
LENGTH WITH SINGLE RING HANDLE GRIP

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