

Chemical Indicator Tape

Multiple widths available for Steam, EO Gas and Hydrogen Peroxide sterilization processes.

Description

Indicator tapes provide immediate identification of processed items while also securing sterilization packs. Color change indicator verifies exposure to sterilization process. Hydrogen peroxide tape comes in its own dispensing box for ease of use. Blue steam tape has been specially designed with aggressive adhesive for disposable CSR wraps. LF2 steam tape is not manufactured with latex and lead. Tape dispenser can hold up to 2 rolls of 1" tape and has a weighted base.

Features & Benefits

- Easy read color change indicator
- 60 yards per roll
- Available in multiple widths
- Individually packaged



ORDERING INFORMATION

Cat. No.	Description	Packaging
HT-048	3/4" Hydrogen Peroxide Indicator Tape	24/case
LF2-048	3/4" LF2 Steam Indicator Tape	24rls/cs
LF2-036	1" LF2 Steam Indicator Tape	18rls/cs
ST-072	1/2" Steam Indicator Tape	36rls/cs
ST-048	3/4" Steam Indicator Tape	24rls/cs
ST-036	1" Steam Indicator Tape	18rls/cs
BT-048	3/4" Blue Steam Indicator Tape	24rls/cs
BT-036	1" Blue Steam Indicator Tape	18rls/cs
GT-048	3/4" EO Gas Indicator Tape	24rls/cs
GT-036	1" EO Gas Indicator Tape	18rls/cs
TD-001	Tape Dispenser	each